



## Real time qPCR protocol - FSQ-301 -

### 【Purpose】

Removal of contaminated genomic DNA in total RNA prior to the cDNA synthesis.

### 【Material】

cDNA synthesis reagent

- ReverTra Ace<sup>®</sup> qPCR RT Master Mix with gDNA Remover (FSQ-301)

Real-time PCR reagent

- THUNDERBIRD<sup>®</sup> SYBR<sup>®</sup> qPCR Mix (QPS-201)

### 【Method】

#### 1. RNA extraction from cells

Total RNA is extracted and purified from HeLa cells using RNA isolation kits

HeLa cell RNA (        ng/uL)

#### 2. cDNA synthesis

cDNA is synthesized using FSQ-301.

Negative control is added No-RT control.


Dilution of RNA (        ng/uL) → 100ng/uL (RNA: 50 uL + Water:        uL)

<b>FSQ-301</b>		
	(+)	(-)
4 × DN Master Mix	4	4
RNA (100ng/uL)	5	5
Water	7	7
<hr/>		
total volume	16	16
↓		
37°C, 5min (4°C, hold)		
↓		
	(+)	(-)
Reaction mix	16	16
5 × RT Master Mix II	4	0
5 × RT Master Mix II no-RT Contr.	0	4
<hr/>		
total volume	20	20
↓		
37°C, 15min 50°C, 5min 98°C, 5min (4°C, hold)		
↓		
cDNA		

### 3. Real-time PCR

RNA detection using QPS-201

QPS-201		
	x1	x15
QPS-201	10	150
50 × ROX	0.4	6
10uM Forward Primer	0.6	9
10uM Reverse Primer	0.6	9
cDNA	2	-
Water	6.4	96
	20	300



1 FSQ-201(+)

2 FSQ-201(-)

3 FSQ-301(+)

4 FSQ-301(-)


↓ Dispense 18uL reaction mix to each tube  
(The test is performed in duplicate.)

↓ Add 2uL each synthesized cDNA


PCR cycle condition  
95°C, 60sec → 95°C, 15sec / 60°C, 45sec (40 cycles) → Melting Curve

#### 【Primer Design】

##### ① TNF (Homo sapiens tumor necrosis facot)



##### ② ACTB (Homo sapiens actin, beta)



##### ③ AMIGO3 (Homo sapiens adhesion molecule with Ig-like domain 3)

