



Frederick Sanger

DNA Sequencing Service Package

IDR 173,000 for 1 Reaction per Sample
IDR 247,000 for 2 Reactions per Sample

Service Includes:

- Gel Electrophoresis for all received samples
- Sample Pre-treatment, which includes:
PCR Clean- Up using Spin column, or Gel Extraction.

DNA Sequencing Service + PLUS

IDR 90,000 for 1 Reaction per Sample
IDR 173,000 for 2 Reactions per Sample

Customer to provide single band un-purified PCR product or agarose gel slice.

For un-purified long PCR product, we guarantee minimum 1,000 bases from a good quality template in a single reaction.

Requirement:

Minimum 30 μL of PCR products in dH_2O per reaction, with 100 $\text{ng}/\mu\text{L}$ concentration.

Minimum 20 μL of primer in dH_2O per reaction, with 10 $\text{pmol}/\mu\text{L}$ concentration.

Single Pass DNA Sequencing, 96-well format

IDR 6,300,000 for 1 plate or 80 - 95 reactions.

High Throughput (HT) DNA Sequencing. Leave well H12 empty for our control reaction, Each plate comprises of maximum 95 samples with normalized OD.

Turnaround Time: 7 working days

*Price is exclude VAT 10%

