


Accommodates Any Sample



Pathogen Inactivation



Inactivates viruses, bacteria,
yeast & protists

DNA/RNA Shield™



Break the Cold Chain Not the bank!




Transport at ambient
temperatures

Streamlined Purification of DNA & RNA



No reagent removal
Compatible with all
purification kits
(including Quick-RNA™,
Quick-DNA™, etc.)
Fully automatable



Ready for all downstream applications

NGS

ARRAYS

(RT)PCR



DNA/RNA Shield™

DNA/RNA Shield™ adalah solusi preservasi sampel yang melindungi DNA dan RNA dari degradasi tanpa perlu penyimpanan beku. Dengan kemampuan menonaktifkan patogen secara instan, produk ini ideal untuk transportasi dan penyimpanan sampel biologis seperti swab, saliva, darah, serta jaringan. Cocok untuk aplikasi penelitian genomik, diagnostik molekuler, dan biobank.

Informasi Produk:

R1100-50	- 50ml
R1100-250	- 250ml
R1200-25	- 25ml (2X Concentrate)
R1200-125	- 125ml (2X Concentrate)

DNA/RNA Shield™ Collection Devices



Worry free solution for sample collection, virus inactivation, and room temp storage of RNA for over a month!

DNA/RNA Shield™ Collection Devices

DESKRIPSI	UKURAN	CAT NO.
DNA/RNA Shield™ Fecal Collection Tube	10 Pack	R1101
DNA/RNA Shield™ SafeCollect™ Swab Collection	1 Kit	R1161
DNA/RNA Shield™ Collection Tube w/ Swab	50 Pack	R1106
DNA/RNA Shield™ SafeCollect™ Saliva Collection	1 Kit	R1211
DNA/RNA Shield™ Blood Collection Tube	50 Pack	R1150

SCAN QR



UNTUK INFORMASI LEBIH LANJUT