

DNA FROM ANY SAMPLE TYPE


Quick-DNA™ Miniprep Plus Kit




Quick-DNA™ Miniprep Plus Kit

- Purify high-quality DNA easily and reliably from any sample source (biological fluids, cultured/monolayer cells, solid tissues, etc.).
- Zymo-Spin™ technology ensures DNA is ultra-pure and ready for all sensitive downstream applications such as qPCR, DNA-sequencing, arrays, and methylation analysis.

Quick-DNA™ Plus Workflow



Quick-DNA™ Plus Performance



The **Quick-DNA™ Miniprep Plus Kit** yields highly concentrated DNA without RNA contamination when compared to the QIAamp® (Qiagen) utilizing porcine whole blood, HeLa cells, and bovine muscle tissue. The resulting DNA from each sample was analyzed in a 1% (w/v) TAE/agarose/EtBr gel.

Quick-DNA™ Plus Ordering Information

Quick-DNA™ Plus (Proteinase K Included)

Any Sample Type - Tissue, Cells, Whole Blood, etc.

Product	Cat. No.	Size	DNA Recovery	Minimum Elution	Input
Quick-DNA™ Microprep Plus Kit	D4074	50 preps	5 µg	10 µl	≤ 10 ⁶ cells ≤ 5 mg tissue
Quick-DNA™ Miniprep Plus Kit	D4068T D4068 D4069	10 preps 50 preps 200 preps	25 µg	35 µl	≤ 5 x 10 ⁶ cells ≤ 25 mg tissue
Quick-DNA™ Midiprep Plus Kit	D4075	25 preps	125 µg	200 µl	≤ 3 x 10 ⁷ cells ≤ 125 mg tissue
Quick-DNA™ 96 Plus Kit	D4070 D4071	2 x 96 preps 4 x 96 preps	5 µg	15 µl	≤ 10 ⁶ cells ≤ 5 mg tissue
Quick-DNA™ Magbead Plus Kit	D4081 D4082	1 x 96 preps 4 x 96 preps	10 µg	75 µl	≤ 3 x 10 ⁶ cells ≤ 25 mg tissue



PT. Genetika Science Indonesia

Green Lake City, Rukan Great Wall Blok C No. 19-21,
Kel. Petir, Kec. Cipondoh, Kota Tangerang,
Banten 15147, Indonesia.

SCAN HERE

