



Genetika
Science



BIOFIL®



JET BIOFIL has consistently **produced high-quality cell culture products** over the more than two decades to ensure agree the best repeatability and the most reliable results in studies, and our products are therefore widely used by researchers. The cell culture products of JET BIOFIL are suitable for most applications and can meet your different demands for cell culture.

CELL & TISSUE CULTURE FLASK

Product Description	Catalog Number	Volume (mL)	Type of Cap	Culture Area (cm ²)	Recommended Working Volume (mL)	Qty. Per Bag
25 mL (volume) TC-Treated	TCF011025	25 mL	Plug Seal	12.5 cm ²	2.5 -3.8 mL	10 Qty

CELL & TISSUE CULTURE DISH

Product Description	Catalog Number	Diameter (mm)	Height (mm)	Culture Area (cm ²)	Recommended Working Volume (mL)	Qty. Per Bag
TC-Treated, 60 mm	TCD010060	60 mm	17.3 mm	21.2 cm ²	4-7 mL	10 Qty

CELL & TISSUE CULTURE DISH

Product Description	Catalog Number	Surface Type	Total of Well	Culture Area (cm ²)	Recommended Working Volume of single Well (mL)	Qty. Per Bag
6-well plates	TCP011006	TC-treated	6 Well	9.6 cm ²	1.9-2.9 mL	1 Qty
24-well plates	TCP011024	TC-treated	24 Well	1.93 cm ²	0.38-0.57 mL	1 Qty
96-well plates	TCP011096	TC-treated	96 Well	0.33 cm ²	0.075-0.2 mL	1 Qty



**BIOFIL®****YOUR RELIABLE PARTNER IN LIFE SCIENCE****SEROLOGICAL PIPETS**

Product Description	Catalog Number	Graduation (mL)	Volume (mL)	Length (mm)	Qty Per Bag
Serological Pipets - Sterile	GSP010010	1/10 mL	10 mL	346.3 mm	1 Qty
Serological Pipets - Sterile	GSP010025	2/10 mL	25 mL	308.5 mm	1 Qty

PIPETTE MICRO TIPS

Product Description	Catalog Number	Type of Tips	Capacity (µL)	Type of Package	Qty Per Rack Box
ZEROTIP® Pipette Micro Tips 0.1-10 µL	PMT251010	Sterile No Filter	0.1-10 µL	Rack Box	96 Qty
ZEROTIP® Pipette Micro Tips 10-200 µL	PMT251200	Sterile No Filter	10-200 µL	Rack Box	96 Qty
ZEROTIP® Pipette Micro Tips 100-1000 µL	PMT251000	Sterile No Filter	100-1000 µL	Rack Box	96 Qty

MICRO CENTRIFUGE TUBE

Product Description	Catalog Number	Capacity (mL)	Bottom	Maximum RCF (xg)	Color	Qty Per Bag
1.5 mL Micro Centrifuge Tubes with Snap Flat Cap (non-sterile)	CFT000015	1.5 mL	Conical	25,000 (xg)	Natural	500 Qty

CENTRIFUGE TUBES

Product Description	Catalog Number	Capacity (mL)	Bottom	Maximum RCF (xg)	Package	Qty Per Bag
15 mL, Sterile, in Resealable Bag	CFT011150	15 mL	Conical	12,000 xg	Re-sealable Bag	25 Qty
50 mL, Sterile, Self Standing, in Resealable Bag	CFT111500	50 mL	Self Standing	6,000 xg	Re-sealable Bag	25 Qty

CRYOGENIC VIALS (NON-STERILE)

Product Description	Catalog Number	Capacity (mL)	Bottom	Cap Color	Color	Qty Per Bag
2 mL Cryogenic Vials with Flat Cap	FCT511020	2.0 mL	Self Standing	Natural	Natural	50 Qty

CELL SCRAPER

Product Description	Catalog Number	Total Length (cm)	Blade (mm)	Type of Package	Qty Per Bag
Cell Scraper 30 cm with 19.5 mm Blade Rotatable	CSC211030	30 cm	19.5 mm	Blades/PE; Handle/ABS	1 Qty

**SCAN QR FOR FUTHER INFORMATIONS & DETAILS**